

# 1-STEP 1-MINUTE

## Safetec SaniZide Pro 1<sup>®</sup> Spray (1-Minute Surface Disinfectant)

Safetec's SaniZide Pro 1<sup>®</sup> is a broad spectrum, EPA Registered surface disinfectant that combines cleaning and disinfecting in one simple step with its ready-to-use, alcohol-based, hospital-grade formula. SaniZide Pro 1<sup>®</sup> makes staying compliant even easier by reducing the risk of cross-contamination with true, one-minute kill claims that are effective against 48 microorganisms and safe for repeated use on hard, non-porous surfaces.



| Item# | Description                                | Case Qty. | Case Wt. |
|-------|--|-----------|----------|
| 35910 | SaniZide Pro 1 <sup>®</sup> Spray (32 oz.) | 6/case    | 14 lbs.  |

EPA Reg. No. 88494-3-67161

**Effective Against MRSA, VRE,  
 Tuberculosis, Influenza, Poliovirus,  
 Norovirus, and SARS-CoV-2 (the  
 virus that causes COVID-19)**



**Medical**



**Dental**



**Veterinary**

### Ideal for use in:

- Hospitals
- Medical or Dental Offices
- Ambulances
- Ambulatory Surgical Centers
- Nursing Homes
- Laboratories
- Pharmacies
- Veterinary Offices
- Other Healthcare Related Institutions

# SaniZide Pro 1® Spray Efficacy Sheet

†Enveloped/lipophilic type virus  
 ††Non-enveloped/hydrophilic type virus

## Test Organism

## Kill Time

|  |          |
|--|----------|
| <i>Acinetobacter baumannii</i> MDR.....  | 1-minute |
| ††Adenovirus.....  | 1-minute |
| †Avian Influenza.....  | 1-minute |
| <i>Bordetella bronchiseptica</i> .....   | 1-minute |
| <i>Bordetella pertussis</i> .....  | 1-minute |
| <i>Burkholderia cepacia</i> .....  | 1-minute |
| <i>Campylobacter jejuni</i> .....  | 1-minute |
| <i>Candida albicans</i> .....  | 1-minute |
| †Canine distemper virus .....  | 1-minute |
| ††Canine parvovirus .....  | 1-minute |
| Carbapenem resistant <i>Enterobacter cloacae</i> (CRE) or NDM-1 Positive .....                             | 1-minute |
| †Coronavirus.....  | 1-minute |
| †Cytomegalovirus .....   | 1-minute |
| <i>Enterobacter aerogenes</i> or <i>Enterobacter cloacae</i> .....   | 1-minute |
| ††Enterovirus EV-D68 .....   | 1-minute |
| <i>Escherichia coli</i> - ESBL Producing MDR or NDM-1 Positive .....                                       | 1-minute |
| <i>Escherichia coli</i> (E. coli) or E. coli O157:H7 .....   | 1-minute |
| †Feline leukemia virus .....   | 1-minute |
| †Hepatitis B Virus (HBV) .....   | 1-minute |
| †Hepatitis C Virus (HCV) .....   | 1-minute |
| †Herpes Simplex Virus.....   | 1-minute |
| †HIV-Aids .....  | 1-minute |
| †Influenza A and/or Pandemic 2009 H1N1 Influenza A, or H3N3 or Hong Kong Flu, or H2N2 .....                | 1-minute |
| †Influenza B .....   | 1-minute |
| <i>Klebsiella pneumoniae</i> .....   | 1-minute |
| <i>Klebsiella pneumoniae</i> -ESBL Producing MDR or KPC-2 or CRKP or NDM-1 Positive.....                   | 1-minute |
| <i>Legionella pneumophila</i> .....  | 1-minute |
| <i>Listeria monocytogenes</i> .....  | 1-minute |
| <i>Mycobacterium bovis</i> (Tuberculosis).....   | 1-minute |
| ††Norovirus .....  | 1-minute |
| ††Poliovirus.....  | 1-minute |
| <i>Proteus vulgaris</i> .....  | 1-minute |
| <i>Pseudomonas aeruginosa</i> .....  | 1-minute |
| †Respiratory Syncytial virus (RSV).....  | 1-minute |
| ††Rhinovirus .....   | 1-minute |
| ††Rotavirus.....   | 1-minute |
| <i>Salmonella enterica</i> .....   | 1-minute |
| <i>Shigella dysenteriae</i> .....  | 1-minute |
| <i>Staphylococcus aureus</i> .....   | 1-minute |
| <i>Staphylococcus aureus</i> - Methicillin resistant MRSA (Healthcare or Community Acquired) or MRSE ..... | 1-minute |
| <i>Staphylococcus aureus</i> - Vancomycin intermediate VISA .....  | 1-minute |
| <i>Staphylococcus aureus</i> - Vancomycin resistant VRSA.....  | 1-minute |
| <i>Streptococcus pneumoniae</i> .....  | 1-minute |
| <i>Streptococcus pyogenes</i> .....  | 1-minute |
| <i>Trichophyton mentagrophytes</i> (Ringworm, Athlete's Foot, Jock Itch) .....                             | 1-minute |
| †Vaccinia virus .....  | 1-minute |
| Vancomycin resistant <i>Enterococcus faecium</i> or <i>faecalis</i> (VRE) .....                            | 1-minute |

The US EPA also recognizes the efficacy of Safetec® SaniZide Pro 1® against the COVID-19 strain of Coronavirus in 1-minute

## May be used on hard, non-porous surfaces in:

**Hospitals and healthcare settings, including:** Ambulances, Anesthesia, CAT Labs, CCU, Central Supply, Clinics, Dental Offices, Dialysis Clinics, Doctor's Offices, Emergency Rooms, Home Health Care, Hospices, Laboratory, Laundry Rooms, Long Term Care Facilities, Newborn Nursery, Nursing Homes, Offices, Operating Rooms, Pediatrics, Physical Therapy, Physician's Offices, Ophthalmic Offices, Radiology, Recovery Rooms, Rehabilitation, Respiratory Therapy, Ultrasound, Transport Vehicles, Waiting Areas, X-Ray. **Critical Care Areas:** Cardiac Care Units, Emergency Rooms, Intensive Care Units, Neonatal Intensive Care Units, Pediatric Intensive Care Units, Surgery, Surgical Prep, and Surgical Recovery Units. **Animal Housing/Holding Areas:** Animal Hospitals, Animal Life Science Laboratories, Animal Transportation Vehicles, Breeding Establishments, Grooming Establishments, Kennels, Pet Shops, Tack Shops, Veterinary Clinics, Veterinary Examination Rooms, Veterinary Operating Rooms, Veterinary Waiting Rooms, Veterinary X-Ray Rooms, Zoos.

SaniZide Pro 1® contains Ethanol (CAS 64-17-5) as a disinfectant, which in California is considered to be carcinogenic to humans when ingested in the form of an alcoholic beverage [IARC 1]. SaniZide Pro 1® also contains the following ingredients which do not appear on any listing of chemicals of concern as designated by the State of California or the European Union: Purified Water (CAS 7732-18-5) as a diluent, Propylene Glycol (CAS 57-55-6) as a solvent, Citric Acid (CAS 77-92-9) as a pH adjuster, Decyl Glucoside (CAS 68515-73-1) as a surfactant, Didecylidimonium Chloride (CAS 7173-51-5) as a disinfectant, Lauryl Glucoside (CAS 110615-47-9) as a surfactant.

# How much surface disinfectant does your facility need?

According to the [AIHA Reopening Guidance Document](#), all contact surfaces should be disinfected regularly, including:

- **Individual workspaces between each shift.**
- **Common area contact surfaces between each use, including but not limited to:**
  - Door handles
  - Elevator buttons
  - Sign-in areas

All calculations are estimates based on kill claims for SaniZide Pro 1® which has a 1-minute contact time.

- **There are 0.0125 oz. in 1 spray of SaniZide Pro 1®**
- **Each SaniZide Pro 1® bottle has 32 oz.**



With this information we can gather:

| Item# | Description                    | Case Qty. |
|-------|--------------------------------|-----------|
| 35910 | SaniZide Pro 1® Spray (32 oz.) | 6/case    |

| Item                              | Number of Sprays | Ounces of SaniZide Pro 1® |
|-----------------------------------|------------------|---------------------------|
| Doorknob                          | 2                | 0.025                     |
| Work Desk                         | 10               | 0.125                     |
| Large conference room table       | 15               | 0.1875                    |
| Light switches                    | 1                | 0.0125                    |
| Elevator buttons                  | 2                | 0.025                     |
| Non-food contact Kitchen surfaces | 10               | 0.125                     |
| Bathroom Surfaces                 | 20               | 0.25                      |
| Sign-in Areas/Desk                | 10               | 0.125                     |

This information is based on the amount of product needed in a year, based on 260 work days a year.

| Surface                     | # Number of Disinfections a Day<br>(Once a Day or After Each Use) | # of Ounces | # of Bottles of<br>SaniZide Pro 1® |
|-----------------------------|---|-------------|------------------------------------|
| Large Conference Room Table | 5   | 244         | 8                                  |
| Sign-in Area                | 5   | 163         | 5                                  |
| 100 Light Switches          | 1   | 325         | 10                                 |
| 100 Doorknobs               | 2   | 1,300       | 40                                 |
| 2 Elevators                 | 100   | 1,250       | 39                                 |
| Unshared Work Desk          | 1   | 32.5        | 1                                  |
| 20 Unshared Work Desks      | 1   | 650         | 20                                 |
| 50 Unshared Work Desks      | 1   | 1,625       | 51                                 |

\*Amount of product needed to keep the surface wet for a full minute varies by climate. Number of sprays are approximate and may vary due to the size of your workspaces and surfaces. The amount of disinfecting necessary will vary by building structure, and daily office activities. Always follow the directions on the label of your surface disinfectant.



# DISINFECTANT SPRAY EFFICACY COMPARISON CHART



**Easy To Use, No Expiration, Hospital-Grade, 1-Min Kill Claims**  
**For use on hard, non-porous, high-contact surfaces**



## May be used on hard, non-porous surfaces in:

**Hospitals and healthcare settings, including:** Ambulances, Anesthesia, CAT Labs, CCU, Central Supply, Clinics, Dental Offices, Dialysis Clinics, Doctor's Offices, Emergency Rooms, Home Health Care, Hospices, Laboratory, Laundry Rooms, Long Term Care Facilities, Newborn Nursery, Nursing Homes, Offices, Operating Rooms, Pediatrics, Physical Therapy, Physician's Offices, Ophthalmic Offices, Radiology, Recovery Rooms, Rehabilitation, Respiratory Therapy, Ultrasound, Transport Vehicles, Waiting Areas, X-Ray. **Critical Care Areas:** Cardiac Care Units, Emergency Rooms, Intensive Care Units, Neonatal Intensive Care Units, Pediatric Intensive Care Units, Surgery, Surgical Prep, and Surgical Recovery Units. **Animal Housing/Holding Areas:** Animal Hospitals, Animal Life Science Laboratories, Animal Transportation Vehicles, Breeding Establishments, Grooming Establishments, Kennels, Pet Shops, Tack Shops, Veterinary Clinics, Veterinary Examination Rooms, Veterinary Operating Rooms, Veterinary Waiting Rooms, Veterinary X-Ray Rooms, Zoos.

| Competitors                        | Safetec              | PDI<br>BE THE DIFFERENCE    | Metrex<br>Protecting People | THE CLOROX COMPANY | Diversey<br>for a cleaner, healthier future | MICRO-SCIENTIFIC |
|------------------------------------|----------------------|-----------------------------|-----------------------------|--------------------|---|------------------|
| Active Ingredient                  | Alcohol-Quat         | Alcohol-Quat                | Alcohol-Quat                | Alcohol-Quat       | Alcohol-Quat                                | Alcohol-Quat     |
| Company Products                   | SaniZide Pro 1 Spray | Sani-Prime Germicidal Spray | CaviCide 1 Spray            | Disinfecting Spray | Virex TB                                    | Opti-Cide3       |
| <b>Enveloped Viruses</b>           |                      |                             |                             |                    |   |                  |
| Avian Influenza                    | <b>1-minute</b>      | 1-minute                    |                             | < 1-minute         | 3-minutes                                   | 2-minutes        |
| Coronavirus                        | <b>1-minute</b>      | 1-minute                    | 1-minute                    | < 1-minute         | 2-minutes                                   |                  |
| Cytomegalovirus                    | <b>1-minute</b>      | 1-minute                    |                             | < 1-minute         |   |                  |
| Hepatitis B Virus                  | <b>1-minute</b>      | 1-minute                    | 1-minute                    | < 1-minute         | 1-minute                                    | 2-minutes        |
| Hepatitis C Virus                  | <b>1-minute</b>      | 1-minute                    | 1-minute                    | < 1-minute         | 1-minute                                    | 2-minutes        |
| Herpes Simplex                     | <b>1-minute</b>      | 1-minute                    | 1-minute                    | < 1-minute         | 3-minutes                                   | 2-minutes        |
| HIV-Aids                           | <b>1-minute</b>      | 1-minute                    | 1-minute                    | < 1-minute         | 1-minute                                    | 2-minutes        |
| Influenza A                        | <b>1-minute</b>      | 1-minute                    | 1-minute                    | < 1-minute         | 3-minutes                                   | 2-minutes        |
| Influenza B                        | <b>1-minute</b>      | 1-minute                    |                             | < 1-minute         |   |                  |
| Respiratory Syncytial virus (RSV)  | <b>1-minute</b>      | 1-minute                    |                             | < 1-minute         | 3-minutes                                   |                  |
| Canine distemper virus             | <b>1-minute</b>      |                             |                             |                    |   | 3-minutes        |
| Feline leukemia virus              | <b>1-minute</b>      |                             |                             |                    |   | 3-minutes        |
| Vaccinia virus                     | <b>1-minute</b>      | 1-minute                    |                             |                    |   |                  |
| <b>Large Non-Enveloped Viruses</b> |                      |                             |                             |                    |   |                  |
| Adenovirus                         | <b>1-minute</b>      | 1-minute                    | 3-minutes                   | 10-minutes         | 3-minutes                                   |                  |
| Rotavirus                          | <b>1-minute</b>      | 1-minute                    | 1-minute                    | 10-minutes         | 3-minutes                                   | 2-minutes        |
| <b>Small Non-Enveloped Viruses</b> |                      |                             |                             |                    |   |                  |
| Norovirus                          | <b>1-minute</b>      |                             | 1-minute                    | 10-minutes         | <1-minute                                   |                  |
| Poliovirus                         | <b>1-minute</b>      |                             |                             | 10-minutes         | 3-minutes                                   |                  |
| Rhinovirus                         | <b>1-minute</b>      | 1-minute                    |                             | 10-minutes         | 3-minutes                                   | 2-minutes        |
| Enterovirus EV-D68                 | <b>1-minute</b>      | 1-minute                    |                             |                    |   |                  |
| Canine parvovirus                  | <b>1-minute</b>      |                             |                             |                    | 3-minutes                                   |                  |
| <b>Mycobacterium</b>               |                      |                             |                             |                    |   |                  |
| Mycobacterium bovis (Tuberculosis) | <b>1-minute</b>      | 1-minute                    | 1-minute                    | 5-minutes          | 5-minutes                                   | 2-minutes        |
| <b>Encapsulated Bacteria</b>       |                      |                             |                             |                    |   |                  |
| Klebsiella pneumoniae              | <b>1-minute</b>      | 1-minute                    | 1-minute                    | 3-minutes          | 3-minutes                                   |                  |
| Streptococcus pneumoniae           | <b>1-minute</b>      | 1-minute                    |                             | 3-minutes          |   |                  |
| Pseudomonas aeruginosa             | <b>1-minute</b>      | 1-minute                    | 1-minute                    | 3-minutes          | 3-minutes                                   | 2-minutes        |
| Bordetella bronchiseptica          | <b>1-minute</b>      |                             |                             |                    |   | 3-minutes        |
| Bordetella pertussis               | <b>1-minute</b>      | 1-minute                    |                             |                    |   |                  |

| Competitors   |  |  |  |  |  |  |
|---|--|--|---|--|--|--|
| <b>Drug Resistant Bacteria</b>  |  |  |   |  |  |  |
| Acinetobacter baumannii MDR   | 1-minute   | 1-minute   | 1-minute  |  |  |  |
| Carbapenum resistant Enterobacter cloacae (CRE) or NDM-1 Positive                             | 1-minute   | 1-minute   |   |  |  |  |
| Escherichia coli - ESBL Producing MDR or NDM-1 Positive                                       | 1-minute   | 1-minute   | 1-minute  | 3-minutes  |  |  |
| Vancomycin resistant Enterococcus faecium or faecalis (VRE)                                   | 1-minute   | 1-minute   | 1-minute  | 3-minutes  | 3-minutes  | 2-minutes  |
| Klebsiella pneumoniae-ESBL Producing MDR or KPC-2 or CRKP or NDM-1 Positive                   | 1-minute   | 1-minute   | 1-minute  | 3-minutes  |  | 2-minutes  |
| Staphylococcus aureus - Methicillin resistant MRSA (Healthcare of Community Acquired) or MRSE | 1-minute   | 1-minute   | 1-minute  | 3-minutes  | 3-minutes  | 2-minutes  |
| Staphylococcus aureus - Vancomycin resistant VRSA   | 1-minute   | 1-minute   |   |  |  |  |
| Staphylococcus aureus - Vancomycin Intermediate VISA  | 1-minute   | 1-minute   | 1-minute  |  |  |  |
| <b>Gram Positive Bacteria</b>   |  |  |   |  |  |  |
| Listeria monocytogenes  | 1-minute   | 1-minute   |   | 3-minutes  | 3-minutes  | 2-minutes  |
| Staphylococcus aureus   | 1-minute   | 1-minute   | 1-minute  | 3-minutes  | 3-minutes  | 2-minutes  |
| Streptococcus pyogenes  | 1-minute   | 1-minute   |   | 3-minutes  | 3-minutes  |  |
| <b>Gram Negative Bacteria</b>   |  |  |   |  |  |  |
| Burkholderia cepacia  | 1-minute   | 1-minute   | 1-minute  | 3-minutes  | 3-minutes  |  |
| Campylobacter jejuni  | 1-minute   | 1-minute   |   | 3-minutes  | 3-minutes  |  |
| Enterobacter aerogenes or Enterobacter cloacae  | 1-minute   | 1-minute   | 1-minute  | 3-minutes  | 3-minutes  | 2-minutes  |
| Escherichia coli (E. coli) or E. coli O157:H7   | 1-minute   | 1-minute   |   | 3-minutes  | 3-minutes  | 2-minutes  |
| Legionella pneumophila  | 1-minute   | 1-minute   |   |  | 3-minutes  |  |
| Proteus vulgaris  | 1-minute   | 1-minute   |   | 3-minutes  | 3-minutes  |  |
| Salmonella enterica   | 1-minute   | 1-minute   | 1-minute  | 3-minutes  | 3-minutes  | 2-minutes  |
| Shigella dysenteriae  | 1-minute   | 1-minute   |   | 3-minutes  | 3-minutes  |  |
| <b>Pathogenic Fungi</b>   |  |  |   |  |  |  |
| Trichophyton mentagrophytes   | 1-minute   | 1-minute   | 1-minute  | 1-minute   | 3-minutes  | 2-minutes  |
| Candida albicans  | 1-minute   | 1-minute   | 1-minute  | 1-minute   | 3-minutes  |  |
| <b>SARS-CoV-2 (COVID-19)</b>  |  |  |   |  |  |  |
| On EPA List N   | ✓  | ✓  | ✓   | ✓  | ✓  | ✓  |
| Contact Time to kill SARS-CoV-2 (COVID-19) according to the <a href="#">EPA List N</a>        | 1-minute   | 3-minutes  | 3-minutes   | 10-minutes   | 2-minutes  | 2-minutes  |
| <b>Product Attributes</b>   |  |  |   |  |  |  |
| Will not harm most surfaces   | ✓  | ✓  | ✓   | ✓  | ✓  | ✓  |
| <1 min or 1 minute for contact time   | ✓  | ✓  |   |  |  |  |
| Fragrance Free  | ✓  | ✓  | ✓   |  |  | ✓  |
| Low Toxicity Rating   | ✓  |  |   |  | ✓  | ✓  |
| No Residue  | ✓  |  | ✓   | ✓  |  | ✓  |
| One-Step Clean and Disinfect  | ✓  | ✓  | ✓   | ✓  | ✓  |  |
| Drug-Resistant Kills  | ✓  | ✓  | ✓   | ✓  | ✓  | ✓  |
| Kills Norovirus and Poliovirus in 1 min or less   | ✓  |  |   |  |  |  |
| Made in USA   | ✓  | ✓  | ✓   | ✓  | ✓  | ✓  |