

1-STEP 1-MINUTE

Safetec SaniZide Pro 1® Wipes (1-Minute Surface Disinfectant)

Safetec's SaniZide Pro 1® is a broad spectrum, EPA Registered surface disinfectant that combines cleaning and disinfecting in one simple step with its ready-to-use, alcohol-based, hospital-grade formula. SaniZide Pro 1® makes staying compliant even easier by reducing the risk of cross-contamination with true, one-minute kill claims that are effective against 48 microorganisms and safe for repeated use on hard, non-porous surfaces.

Item#	Description	Case Qty.	Case Wt.
35926	SaniZide Pro 1® Wipes (50 ct. box)*	6/case	12 lbs.
35923	SaniZide Pro 1® Wipes (150 ct. canister)**	12/case	23 lbs.

EPA Reg. No. 88494-4-67161

*Wipe size: 7.5" x 8"

**Wipe size: 6" x 6.75"

Ideal for use in:

- Hospitals
- Medical or Dental Offices
- Ambulances
- Ambulatory Surgical Centers
- Nursing Homes
- Laboratories
- Pharmacies
- Other Healthcare Related Institutions



**Effective Against MRSA, VRE,
 Tuberculosis, Influenza, Poliovirus,
 Norovirus, and SARS-CoV-2 (the
 virus that causes COVID-19)**



Medical



Dental



Made in the USA

SaniZide Pro 1® Wipes Efficacy Sheet

†Enveloped/lipophilic type virus
 ††Non-enveloped/hydrophilic type virus

Test Organism

Kill Time

<i>Acinetobacter baumannii</i> MDR.....	1-minute
††Adenovirus.....	1-minute
†Avian Influenza.....	1-minute
<i>Bordetella bronchiseptica</i>	1-minute
<i>Bordetella pertussis</i>	1-minute
<i>Burkholderia cepacia</i>	1-minute
<i>Campylobacter jejuni</i>	1-minute
<i>Candida albicans</i>	1-minute
†Canine distemper virus	1-minute
††Canine parvovirus	1-minute
Carbapenem resistant <i>Enterobacter cloacae</i> (CRE) or NDM-1 Positive	1-minute
†Coronavirus.....	1-minute
†Cytomegalovirus	1-minute
<i>Enterobacter aerogenes</i> or <i>Enterobacter cloacae</i>	1-minute
††Enterovirus EV-D68	1-minute
<i>Escherichia coli</i> - ESBL Producing MDR or NDM-1 Positive	1-minute
<i>Escherichia coli</i> (E. coli) or E. coli O157:H7	1-minute
†Feline leukemia virus	1-minute
†Hepatitis B Virus (HBV)	1-minute
†Hepatitis C Virus (HCV)	1-minute
†Herpes Simplex Virus.....	1-minute
†HIV-Aids	1-minute
†Influenza A and/or Pandemic 2009 H1N1 Influenza A, or H3N3 or Hong Kong Flu, or H2N2	1-minute
†Influenza B	1-minute
<i>Klebsiella pneumoniae</i>	1-minute
<i>Klebsiella pneumoniae</i> -ESBL Producing MDR or KPC-2 or CRKP or NDM-1 Positive.....	1-minute
<i>Legionella pneumophila</i>	1-minute
<i>Listeria monocytogenes</i>	1-minute
<i>Mycobacterium bovis</i> (Tuberculosis).....	1-minute
††Norovirus	1-minute
††Poliovirus.....	1-minute
<i>Proteus vulgaris</i>	1-minute
<i>Pseudomonas aeruginosa</i>	1-minute
†Respiratory Syncytial virus (RSV).....	1-minute
††Rhinovirus	1-minute
††Rotavirus.....	1-minute
<i>Salmonella enterica</i>	1-minute
<i>Shigella dysenteriae</i>	1-minute
<i>Staphylococcus aureus</i>	1-minute
<i>Staphylococcus aureus</i> - Methicillin resistant MRSA (Healthcare or Community Acquired) or MRSE	1-minute
<i>Staphylococcus aureus</i> - Vancomycin intermediate VISA	1-minute
<i>Staphylococcus aureus</i> - Vancomycin resistant VRSA.....	1-minute
<i>Streptococcus pneumoniae</i>	1-minute
<i>Streptococcus pyogenes</i>	1-minute
<i>Trichophyton mentagrophytes</i> (Ringworm, Athlete's Foot, Jock Itch)	1-minute
†Vaccinia virus	1-minute
Vancomycin resistant <i>Enterococcus faecium</i> or <i>faecalis</i> (VRE)	1-minute

The US EPA also recognizes the efficacy of Safetec® SaniZide Pro 1® against the COVID-19 strain of Coronavirus in 1-minute

May be used on hard, non-porous surfaces in:

Hospitals and healthcare settings, including: Ambulances, Anesthesia, CAT Labs, CCU, Central Supply, Clinics, Dental Offices, Dialysis Clinics, Doctor's Offices, Emergency Rooms, Home Health Care, Hospices, Laboratory, Laundry Rooms, Long Term Care Facilities, Newborn Nursery, Nursing Homes, Offices, Operating Rooms, Pediatrics, Physical Therapy, Physician's Offices, Ophthalmic Offices, Radiology, Recovery Rooms, Rehabilitation, Respiratory Therapy, Ultrasound, Transport Vehicles, Waiting Areas, X-Ray. **Critical Care Areas:** Cardiac Care Units, Emergency Rooms, Intensive Care Units, Neonatal Intensive Care Units, Pediatric Intensive Care Units, Surgery, Surgical Prep, and Surgical Recovery Units.

SaniZide Pro 1® contains Ethanol (CAS 64-17-5) as a disinfectant, which in California is considered to be carcinogenic to humans when ingested in the form of an alcoholic beverage [IARC 1]. SaniZide Pro 1® also contains the following ingredients which do not appear on any listing of chemicals of concern as designated by the State of California or the European Union: Purified Water (CAS 7732-18-5) as a diluent, Propylene Glycol (CAS 57-55-6) as a solvent, Citric Acid (CAS 77-92-9) as a pH adjuster, Decyl Glucoside (CAS 68515-73-1) as a surfactant, Didecyldimmonium Chloride (CAS 7173-51-5) as a disinfectant, Lauryl Glucoside (CAS 110615-47-9) as a surfactant.

DISINFECTANT WIPE EFFICACY COMPARISON CHART



**Easy To Use, No Expiration, Hospital-Grade,
1-Min Kill Claims**

*(Suggested use for disinfecting: ambulance equipment,
patient care surfaces and medical devices)*



Competitors	Safetec	PDI BE THE DIFFERENCE	Metrex Protecting People	THE CLOROX COMPANY	PALMERO HEALTHCARE Protecting Clinicians, Patients, and the Practice	MICRO-SCIENTIFIC
Active Ingredient	Alcohol-Quat	Alcohol-Quat	Alcohol-Quat	Alcohol-Quat	Alcohol-Quat	Alcohol-Quat
Company Products	SaniZide Pro 1 Wipes	Sani-Cloth Prime	CaviWipe 1	Multi-Surface Quat Alcohol Cleaner Disinfectant Wipes	Discide Ultra	Opti-Cide3
Enveloped Viruses						
Avian Influenza	1-minute	1-minute		2-minute	< 1-minute	2-minute
Coronavirus	1-minute	1-minute	1-minute		< 1-minute	
Cytomegalovirus	1-minute	1-minute				
Hepatitis B Virus	1-minute	1-minute	1-minute	2-minute	1-minute	2-minute
Hepatitis C Virus	1-minute	1-minute	1-minute	2-minute	1-minute	2-minute
Herpes Simplex	1-minute	1-minute	1-minute	2-minute	1-minute	2-minute
HIV-Aids	1-minute	1-minute	1-minute	2-minute	1-minute	2-minute
Influenza A	1-minute	1-minute	1-minute	2-minute	< 1-minute	2-minute
Influenza B	1-minute	1-minute				
Respiratory Syncytial virus (RSV)	1-minute	1-minute			1-minute	
Canine distemper virus	1-minute			3-minute		3-minute
Feline leukemia virus	1-minute			3-minute		3-minute
Vaccinia virus	1-minute	1-minute			1-minute	
Large Non-Enveloped Viruses						
Adenovirus	1-minute	1-minute	3-minutes		1-minute	
Rotavirus	1-minute	1-minute		2-minute		2-minute
Small Non-Enveloped Viruses						
Norovirus	1-minute					
Poliovirus	1-minute					
Rhinovirus	1-minute	1-minute		2-minute		2-minute
Enterovirus EV-D68	1-minute	1-minute				
Canine parvovirus	1-minute			5-minute		5-minute
Mycobacterium						
Mycobacterium bovis (Tuberculosis)	1-minute	1-minute	1-minute	2-minute	1-minute	2-minute
Encapsulated Bacteria						
Klebsiella pneumoniae	1-minute	1-minute	1-minute			
Streptococcus pneumoniae	1-minute	1-minute				
Pseudomonas aeruginosa	1-minute	1-minute	1-minute	2-minute	1-minute	2-minute
Bordetella bronchiseptica	1-minute			3-minute		3-minute
Bordetella pertussis	1-minute	1-minute				

Competitors	 Safetec	 PDI BE THE DIFFERENCE	 Metrex Protecting People	 THE CLOROX COMPANY	 PALMERO HEALTHCARE Protecting Clinicians, Patients, and the Practice	 MICRO-SCIENTIFIC
Drug Resistant Bacteria						
Acinetobacter baumannii MDR	1-minute	1-minute	1-minute			
Carbapenum resistant Enterobacter cloacae (CRE) or NDM-1 Positive	1-minute	1-minute				
Escherichia coli - ESBL Producing MDR or NDM-1 Positive	1-minute	1-minute	1-minute	2-minute		2-minute
Vancomycin resistant Enterococcus faecium or faecalis (VRE)	1-minute	1-minute	1-minute	2-minute	1-minute	2-minute
Klebsiella pneumoniae-ESBL Producing MDR or KPC-2 or CRKP or NDM-1 Positive	1-minute	1-minute	1-minute	2-minute		2-minute
Staphylococcus aureus - Methicillin resistant MRSA (Healthcare of Community Acquired) or MRSE	1-minute	1-minute	1-minute	2-minute	1-minute	2-minute
Staphylococcus aureus - Vancomycin resistant VRSA	1-minute	1-minute				
Staphylococcus aureus - Vancomycin Intermediate VISA	1-minute	1-minute	1-minute			
Gram Positive Bacteria						
Listeria monocytogenes	1-minute	1-minute		2-minute		2-minute
Staphylococcus aureus	1-minute	1-minute	1-minute	2-minute	1-minute	2-minute
Streptococcus pyogenes	1-minute	1-minute				
Gram Negative Bacteria						
Burkholderia cepacia	1-minute	1-minute	1-minute		1-minute	
Campylobacter jejuni	1-minute	1-minute			1-minute	
Enterobacter aerogenes or Enterobacter cloacae	1-minute	1-minute	1-minute	2-minute		
Escherichia coli (E. coli) or E. coli O157:H7	1-minute	1-minute			1-minute	2-minute
Legionella pneumophila	1-minute	1-minute				
Proteus vulgaris	1-minute	1-minute				
Salmonella enterica	1-minute	1-minute	1-minute	2-minute	1-minute	2-minute
Shigella dysenteriae	1-minute	1-minute		3-minute	3-minute	
Pathogenic Fungi						
Trichophyton mentagrophytes	1-minute	1-minute	3-minutes	2-minute		2-minute
Candida albicans	1-minute	1-minute	1-minute	1-minute	1-minute	
Product Attributes						
Will not harm most surfaces	✓	✓	✓	✓	✓	✓
<1 min or 1 minute for contact time	✓	✓			✓	
Fragrance Free	✓	✓	✓	✓		✓
Low Toxicity Rating	✓			✓		✓
No Residue	✓		✓	✓	✓	✓
One-Step Clean and Disinfect	✓	✓	✓		✓	
Drug-Resistant Kills	✓	✓	✓	✓	✓	✓
Kills Norovirus and Poliovirus in 1 min or less	✓					
Made in USA	✓	✓	✓	✓	✓	✓
Pricing						
Average Price per Wipe	\$0.21/wipe	\$0.25/wipe	\$0.26/wipe	\$0.07/wipe	\$0.16/wipe	\$0.16/wipe